

Endo-Model Knee System

with Segmental Bone Replacement Components and MIRETO Instrument Set

Surgical Technique



Explanat	Explanation of Pictograms									
****	Manufacturer	REF	Article number							
MAT	Material number	RX only	Caution: Federal law restricts this device to sale by or on the order of a physician.							



Endo-Model Knee System

with Segmental Bone Replacement Components and MIRETO Instrument Set

System Description

- 02 Endo-Model Knee Prosthesis System
- 03 LINK PorEx Surface Modification
- 06 Endo-Model Modular Rotating Hinge Knee Prosthesis, Assembly: Modular Stems
- 07 Endo-Model Pure (Non-Rotating) Hinge Knee Prosthesis
- 08 Indications / Contraindications
- 09 MIRETO Instrument Set

Implants

- 10 Endo-Model Rotating Hinge Knee Prosthesis with anti-luxation device
- 11 Endo-Model Rotating Hinge Knee Prosthesis with anti-luxation device and LINK PorEx
- 12 Endo-Model Pure (Non-Rotating) Hinge Knee Prosthesis with hinge axis

Endo-Model Modular Knee Prosthesis System

- 13 Measurements: Joint Components
- Joint Components for Rotating Hinge Version CoCrMo
- 15 Joint Components for Rotating Hinge Version CoCrMo/LINK PorEx
- 16 Joint Components for Pure (Non-Rotating) Hinge Version CoCrMo

Endo-Model Modular Joint Components

- 17 Condylar Replacement
- 18 Intracondylar Version
- 19 Modular Stems, cemented
- 20 Modular Stems, uncemented, conical
- 21 Modular Stems, uncemented, cylindrical
- 22 Femoral Segments Tilastan for Modular Rotating and Pure (Non-Rotating) Hinge Knee Version
- 23 Femoral Segments UHMWPE for Modular Rotating and Pure (Non-Rotating) Hinge Knee Version
- Proximal Tibial Spacers and Segments for Modular Rotating and Pure (Non-Rotating) Hinge Knee Version
- 24+25 UHMWPE
- 26+27 Tilastan
 - Endo-Model Proximal Tibial Spacers for Modular Rotating and Pure (Non-Rotating) Hinge Knee Version
- 28 UHMWPE
- 29+30 Tilastan

31 Centralizers, Patellar Components

- 32 Replacement Sets for Modular Rotating and Pure (Non-Rotating) Hinge Knee Prosthesis
- 33 MIRETO Instrument Set for Endo-Model Knee Prosthesis Systems, case 1 7
- 34 Case 1 General Instruments
- 35 Case 2 Femoral Instruments (2 trays)
- 37 Case 3 Tibial Instruments
- 38 Case 4 Tapered Reamers, conical & cylindrical
- 39 Case 5 Trial Prostheses (2 trays)
- 41 Case 6 Trial Prostheses, cylindrical
- 42 Case 7 Tapered Reamers, conical
- 43 Additional Instrument Set for Endo-Model Rotating Hinge and Modular Knee Prostheses

MIRETO Surgical Technique

- 44 Tibial Preparation
- 45 Tibial Resection
- 53 Femoral Preparation
- 65 Trial Reduction
- 67 Implantation
- 70 Literature

Important Information



Endo-Model Knee Prosthesis Systems

Adherence to the low-friction principle means that the physiological range of motion of this prosthesis is optimally designed thanks to a mounted pivot within the physiological area. The **Endo-Model Rotating Hinge Knee Prosthesis** enables flexion of the joint up to 142°. The joint kinematics also comprise physiological rotation which, given the special shape of the tibial contact surface, ensure elastic force transmission. The hinge knee prosthesis only permits flexion of the joint up to 142° without rotation.

With every step, and especially in the case of falls, torsional stresses are transmitted to the prosthesis anchorage which impact negatively upon the lifespan of the anchorage in the long term. Constructively generated, elastic force transmission ensures preservation at the prosthesis/bone cement and bone cement/bone interfaces. Because of the favorable dimensions of the **Endo-Model Rotating Hinge Knee Prosthesis**, only minimal resection of 14 mm is necessary on the tibia-femur joint plane. The medium-sized intracondylar component is only 30 mm wide. As a rule, there is thus less resection than with a total knee prosthesis. This is a key, positive factor with respect to subsequent revision surgery. The ideal design and dimensions of the knee prosthesis offer a good overview of the operating field. The femoral and tibial components are simply pushed together and the UHMWPE tibial plateau is positioned using a special instrument. Both components are linked by the plateau so as to prevent luxation and without reducing the motional and rotational sequences. The **Pure (Non-Rotating) Hinge Knee Prosthesis** is linked by the axis mechanism.

Flexion and rotational movement of the rotating hinge knee prosthesis is achieved by means of a cross joint. Hyperextension amounts to 2°. The compromise axis lies in the region of the physiological pivot point. Flexion up to 142° is possible. With arthroplasty knee replacement, advancement of the patella or of the patellar bearing surface is often seen. By moving the femoral component posteriorly relative to the tibial axis, physiological movement is also retained for the femoropatellar joint. This protects against progression of retropatellar arthrosis. Rotation of the prosthesis ends in extension by form closure, which ensures a secure standing position. Rotation increases continuously with greater flexion. This rotation is limited primarily by the capsular ligament apparatus. The shape of



the surfaces which are in contact with each other means that further rotation is cushioned elastically by the body weight's bearing down on the joint. The femoral component of the **Endo-Model Rotating Hinge Knee Prosthesis** has a physiological valgus angle of 6°.

Both prosthesis components are broadly supported on their corresponding joint surfaces, such that the compressive strength of the cancellous bone vis-à-vis the femur and tibia is not exceeded. The runnerinspired form of the femoral component is based on the anatomical conditions. The ventral depression provides a smooth transition to the bony patellar bearing surface. The modular prosthesis stems are available for the Endo-Model Modular Rotating Hinge Knee Prosthesis both as a cemented version, without structuring, and with longitudinal structuring for uncemented implantation. To achieve a central position within the medullary cavity, the tips of the cemented stems are fitted with star shaped UHMWPE caps. Direct contact of the metal stems with the inner wall of the bone is thereby prevented. The stems are supplied in lengths of 50 mm up to 280 mm. Special femoral and tibial segments and spacers made of UHMWPE and Tilastan are available for revision surgery of total knee prostheses, to reconstruct condyles and the joint line as well as for tumor cases (resections). It must be noted here that these segments may only be used in combination with corresponding longer stems.

LINK PorEx Surface Modification

LINK PorEx (TiNbN = Titanium Niobium Nitride) Surface Modification

The LINK PorEx surface modification results in a ceramic-like surface, which significantly reduces chrome and nickel ions release.¹

Thanks to its outstanding hardness, abrasion properties similar to ceramics and larger wetting angle, in contact with liquids the LINK PorEx surface against UHMWPE has a lower coefficient of friction compared to CoCrMo surfaces.¹


¹ Internal technical report: Study of the influence of TiNbN-coating on the ion release of CoCrMo-alloys in SBF buffer simulator testing







Endo-Model Knee Prosthesis Systems

Adherence to the low-friction principle means that the physiological range of movement of the Rotating Hinge Knee Prosthesis is optimally designed thanks to a mounted pivot within the physiological area. Flexion and rotational movement of the Rotating Hinge Knee Prosthesis is achieved by means of a cross joint.








The extent of free rotation as a function of flexion and the constructively generated smoothly slowed down rotation area are shown as the hatched area.

Engelbrecht, E.: Die Rotationsendoprothese des Kniegelenks (Rotation prosthesis for the knee joint), Springer-Verlag 1984, ISBN: 978-3-642-69819-4 (print), 978-3-642-69818-7 (online)






Endo-Model Knee Prosthesis Systems


Assembly: Plateau with anti-luxation device




After cementation of tibial and femoral components, the UHMWPE plateau is removed from the tibial tray by loosening the trial screw. With the knee in flexion, the upper and lower components are assembled.



The tibial plateau is attached to the inserter and slid between the femoral and tibial components into the joint so that the plateau chamber grips over the flange. It must be ensured that the dovetail-shaped incision (fig. **A**) on the bottom of the UHMWPE plateau locks into the peripheral groove on the metal tibial support.







The UHMWPE plateau is pressed down and fixed into place by the self-locking fixation screw.

Implanted Endo-Model Modular Knee Prosthesis.





Endo-Model Modular Rotating Hinge Knee Prosthesis

Assembly: Modular Stems

The modular stems are joined by a cone assembly. To ensure rotational stability, the stem has two opposing flanges which are inserted in the medial and lateral grooves on the femoral / tibial components.


The current version features a 6 mm groove for attachment of modular stems with female taper and flanges of 3 mm or 6 mm. When attaching modular stems with 3 mm flanges, the stem must be orientated on the taper so that the threaded hole for the counter screw is not obscured (**A**). To this end, the alignment aid (15-6096/00) for modular stems is used (**B**). Modular stems with 6 mm flanges cannot be combined with Endo-Model implants with 3 mm grooves.





The cone-shaped tip of the locking screw (2) in the taper (3) of the tibial or femoral component pushes the stem firmly against the taper as it is screwed further in (1). A counter screw (4) prevents the locking screw from becoming loose. The screw is secured from medial.

3





Endo-Model Pure (Non-Rotating) Hinge Knee Prosthesis

The external shape, dimensions and sizes of the **Endo-Model Pure (Non-Rotating) Hinge Knee Prosthesis** correspond to those of the Endo-Model Rotating Hinge Knee Prosthesis. As the implant beds required for the Pure and Rotating Hinge versions are identical, the decision as to whether to use a Rotating Hinge or a more stabilizing Pure Hinge Knee Prosthesis can be made intraoperatively.

The connecting piece **A** which is securely attached to the Pure Hinge Knee tibial component and links it to the femoral component is drilled through to accommodate the joint axis **B**. The ventral drilled hole **C** is provided for the grub screw **D** whose tip fits into the recess **E** on the axis and locks the latter once the upper and lower components have been joined.

From inside the intracondylar box of the femoral component, polyethylene bearings **F** for the prosthesis axis are pushed into medial and lateral boreholes. The upper and lower prosthesis components are joined by introducing the tibial connecting piece into the intracondylar box of the femoral component, such that the prosthesis axis can be inserted (always from the medial) using the threaded rod. Articulation takes place between the prosthesis axis and the two bearings.

The Endo-Model Pure (Non-Rotating) Hinge Knee Prosthesis is delivered ready assembled in a sterile condition without centering stars. To disassemble, the grub screw **D** is turned counterclockwise. The threaded rod is screwed onto the prosthesis axis **B** which is then pulled out. The bearings **F** of the upper prosthesis component are pushed inward and removed (when the bearings are subsequently refitted, it must be ensured that the open bearing is positioned medially!).

The package contains two sterile trial bearings (not autoclavable). These are inserted into the upper prosthesis component during surgery; after the trial run, they are exchanged for the definitive bearings. These too can be exchanged if necessary in a second intervention.



Specified Indications and Contraindications for Endo-Model Rotating Hinge and Pure (Non-Rotating) Hinge Knee Prostheses	Rotational version	Hinged version	Implant components with TINbN surface modification
General Indications			
• Severe joint diseases with limitation of mobility due to degenerative, rheumatoid or post-traumatic arthrosis or arthritis. Joint fractures which disallow an osteosynthetic reconstruction.	Х	Х	Х
Indications			
Bone necroses	Х	Х	Х
Bicondylar arthrosis by partly damaged collateral ligaments	Х	_	Х*
 Bicondylar arthrosis by completely damaged ligaments and muscular instability 	_	Х	X*
Revision after primary total knee replacement	Х	Х	Х
Revision surgery after hinge knee or rotational knee joint	Х	Х	Х
 Revision surgery by insufficient / inadequate bone mass 	Х	Х	Х
Differential Indications			
Arthrosis of patella flange	Х	Х	Х
Valgus/Varus deformities <10°	Х	Х	Х
 Valgus/Varus deformities 10 – 15° 	Х	Х	Х
• Valgus/Varus deformities 15 – 20°	Х	Х	Х
• Valgus/Varus deformities 20 – 30°	_	Х	Х*
Contraindications			
 Acute or chronic infections, local and systemic 	Х	Х	Х
Allergies to (implant) materials	Х	Х	Х
Distinctive muscular, nerve, vascular or other diseases which put the affected limb at risk	Х	Х	Х
 Insufficient/inadequate bone mass or -quality which prevents a stable anchoring of the prosthesis. 	Х	Х	Х
Relative Contraindications			
Adiposity	Х	Х	Х
Insufficient musculature	Х	—	X*
Lacking or foreseeable not assured compliance	Х	Х	Х
Foreseeable overload of joint prosthesis	Х	Х	Х

* dependent on the implant variant

This device is intended for cemented use only unless a uncemented modular stem is indicated for use.

These indications/contraindications refer to standard cases. The ultimate decision on whether or not an implant is suitable for a patient must be made by the surgeon based on his/her individual analysis and his/her experience.



MIRETO Instrument-Set

The new **MIRETO Instrument Set** allows users to perform implantation of the Endo-Model Knee Prosthesis System in a safe, bone-conserving, reproducible and precise manner.

Advantages:

All intracondylar Endo-Model standard and modular implants for primary and revision indications can be implanted using the **MIRETO Instrument Set**.

- Minimal number of instrument trays combined with enhanced application modularity.
- Low investment and preparation costs.

MIRETO Instruments were developed with a firm focus on ensuring optimal handling: from storage configuration through to hygienic preparation.


- Accelerated surgical sequence.
- Shorter learning curve.
- Simple hygienic preparation.

The **MIRETO Instrument Set** is characterized by its simple handling, quick and easy assembly/disassembly as well as completely guided instrumentation.

• Reproducible, reliable and precise bone resection by the user.













55

x-small (60) small

65

mediun

(75

large

28 small

30

medium



right

left

Endo-Model Rotating Hinge Knee Prosthesis, with anti-luxation device

MAT EndoDur (CoCrMo), EndoDur – S (CoCrMo), UHMWPE

REF	Side	Size	(R) Radius*
15-8020/11	right	x-small	17 mm
15-8020/12	left	x-small	17 mm
15-8022/11	right	small	20 mm
15-8022/12	left	small	20 mm
15-8024/11	right	medium	23 mm
15-8024/12	left	medium	23 mm
15-8030/11	right	large	25 mm
15-8030/12	left	large	25 mm

* (R) Radius in the sagittal plane: Measured from the center of axis.

Notice: Length specification incl. centralizers* 160 x-small 182 small (185 medium (190 large 6° Valgus Positioning (14) (14) 120 x-small (160) small medium large

* Centralizers are not included in prosthesis packing



Endo-Model Rotating Hinge Knee Prosthesis, with anti-luxation device and LINK PorEx* surface modification

MAT EndoDur (CoCrMo), EndoDur – S (CoCrMo)/LINK PorEx*, UHMWPE

REF	Side	Size	(R) Radius**
15-9020/11	right	x-small	17 mm
15-9020/12	left	x-small	17 mm
15-9022/11	right	small	20 mm
15-9022/12	left	small	20 mm
15-9024/11	right	medium	23 mm
15-9024/12	eft	medium	23 mm
15-9030/11	right	large	25 mm
15-9030/12	left	large	25 mm

** (<u>R) Radius in the sagittal plane</u>: Measured from the center of axis.

Same dimensions as models with anti-luxation device, see pages 10 and 11.



* LINK PorEx: TiNbN = Titanium Niobium Nitride; surface modification (gold color).



Endo-Model Pure (Non-Rotating) Hinge Knee Prosthesis, with hinge axis

MAT EndoDur (CoCrMo), EndoDur- S (CoCrMo), UHMWPE


REF	Side	Size	(R) Radius*
15-2459/11	right	x-small	17 mm
15-2459/12	left	x-small	17 mm
15-2460/11	right	small	20 mm
15-2460/12	left	small	20 mm
15-2461/11	right	medium	23 mm
15-2461/12	eft	medium	23 mm
15-2462/11	right	large	25 mm
15-2462/12	left	large	25 mm

* (R) Radius in the sagittal plane: Measured from the center of axis.




* Centralizers are not included in prosthesis packing








Measurements: Joint Components





Size	A	B	C	FK	X	TK	R
Version	mm	mm	mm	mm	mm	mm	mm
x-small/right	55	42	28	39	50	22	17
x-small/left	55	42	28	39	50	22	17
small/right	60	45	30	42	57	22	20
small/left	60	45	30	42	57	22	20
medium/right	65	45	30	46	62	22	23
medium/left	65	45	30	46	62	22	23
large/right	75	48	35	50	65	22	25
large/left	75	48	35	50	65	22	25



Joint Components Rotating Hinge Knee Version - with anti-luxation device

Joint Components

					0			
MAT EndoDur UHMWPE, Tila		ndoDur – S (C	oCrMo),	MAT EndoDur (CoCrMo), EndoDur – S (CoCrMo), UHMWPE, Tilastan-S				
Modular Join consisting of:	t Componen	t Units		Femoral Corr	ponents:	Tibial Comp	onents:	
REF	Size	Version	Width mm	REF	Version	REF	Version	
15-2815/11	x-small	right	55	15-2810/11	right	15 001 1/01		
15-2815/12	x-small	left	55	15-2810/12	left	15-2814/01	neutral	
15-2816/11	small	right	60	15-2811/11	right	15 0014/00		
15-2816/12	small	left	60	15-2811/12	left	15-2814/02	neutral	
15-2817/11	medium	right	65	15-2812/11	right	15 0014/00	a sudual	
15-2817/12	medium	left	65	15-2812/12	left	15-2814/03	neutral	
15-2818/11	large	right	75	15-2813/11	right	15 0014/04	poutrol	
15-2818/12	large	left	75	15-2813/12	left	15-2814/04	neutral	


Screws to secure the taper assembly between joint component and stem:

A pointed stem locking screw is already located inside the taper of each joint component. The inside packing unit of each joint component includes a counter screw (+ replacement screw) to secure the stem locking screw.



Joint Components Rotating Hinge Knee Version – with LINK PorEx* and anti-luxation device

Joint Components with LINK PorEx*



Screws to secure the taper assembly between joint component and stem:

left

75

large

A pointed stem locking screw is already located inside the taper of each joint component. The inside packing unit of each joint component includes a counter screw (+ replacement screw) to secure the stem locking screw.

15-3818/12



Joint Components Pure (Non-Rotating) Hinge Knee Version

Joint Components

	Fa	0						
MAT EndoDur UHMWPE, Tila		ndoDur – S (C	oCrMo),	MAT EndoDur (CoCrMo), EndoDur – S (CoCrMo), UHMWPE, Tilastan-S				
Modular Join consisting of:	t Componen	t Units	-	Femoral Com	ponents:	Tibial Comp	onents:	
REF	Size	Version	Width mm	REF	Version	REF	Version	
15-2835/11	x-small	right	55	15-2830/11	right	45.0004/04		
15-2835/12	x-small	left	55	15-2830/12	left	15-2834/01	neutral	
15-2836/11	small	right	60	15-2831/11	right	15-2834/02	poutral	
15-2836/12	small	left	60	15-2831/12	left	10-2004/02	neutral	
15-2837/11	medium	right	65	15-2832/11	right	15-2834/03	neutral	
15-2837/12	medium	left	65	15-2832/12	left	10-2004/00	neutrai	
15-2838/11	large	right	75	15-2833/11	right	15-2834/04	neutral	
15-2838/12	large	left	75	15-2833/12	left	10-2004/04	Πουτιαί	

Screws to secure the taper assembly between Joint Component and Stem:

A pointed stem locking screw is already located inside the taper of each joint component. The inside packing unit of each joint component includes a counter screw (+ replacement screw) to secure the stem locking screw.




Endo-Model Modular Joint Components

Condylar Replacement



MAT EndoDur (CoCrMo), EndoDur – S (CoCrMo), UHMWPE, Tilastan-S

Modular Joint Component Units consisting of:			Femoral Compo condylar replace	,	Tibial Components:		
REF	Size	Side	REF	Width mm	REF	Width mm	
15-8521/05	small	right	15-8521/06	60	15,0014/00		
15-8521/07	small	left	15-8521/08	60	15-2814/02	60	
15-8521/09	medium	right	15-8521/10	65	15 0014/00	25	
15-8521/11	medium	left	15-8521/12	65	15-2814/03	65	
15-8521/13	large	right	15-8521/14	75	15 0014/04	75	
15-8521/15	large	left	15-8521/16	75	15-2814/04	75	





Endo-Model Modular Joint Components

Intracondylar Version



MAT CoCrMo, UHMWPE, Tilastan-S

			Femoral Compo intracondylar:	nents,	Tibial Components:		
REF	Size	Side	REF	Width mm	REF	Width mm	
15-8521/25	small	right	15-8521/26	60	15-2814/02	60	
15-8521/27	small	left	15-8521/28	60	15-2014/02		
15-8521/29	medium	right	15-8521/30	65	15,0014/00	0E	
15-8521/31	medium	left	15-8521/32	65	15-2814/03	65	
15-8521/33	large	right	15-8521/34	75	15 0014/04	75	
15-8521/35	large	left	15-8521/36	75	15-2814/04	75	


The joint components are equipped with an anti-luxation device. Femoral components feature female tapers. The femoral components with female taper are not compatible with the femoral segments size 1 to 4!

Femoral component: Width Tibial component: Width



Endo-Model Modular Stems, cemented






Modular Stems, cemented

MAT EndoDur – S (CoCrMo), EndoDur – S (CoCrMo)/LINK PorEx*

			Assembly length**					
REF	REF		Tibia	Femur				
EndoDur – S (CoCrMo)	EndoDur – S (CoCrMo)/ LINK PorEx [*]	L mm	LT mm	LF1 mm	LF2 mm	LF3 mm	LF4 mm	
15-2950/01	15-3950/01	50	87	104	107	111	114	
15-2950/02	15-3950/02	80	117	134	137	141	144	
15-2950/03	15-3950/03	95	132	149	152	156	159	
15-2950/04	15-3950/04	120	157	174	177	181	184	
15-2950/05	15-3950/05	135	172	189	192	196	199	
15-2950/06	15-3950/06	160	197	214	217	221	224	
15-2950/07	15-3950/07	200	237	254	257	261	264	
15-2950/08	15-3950/08	240	277	294	297	301	304	
15-2950/09	15-3950/09	280	317	334	337	341	344	



* LINK PorEx: TiNbN = Titanium Niobium Nitride; surface modification (gold color). Available on request.

** Assembly length incl. centering star unit joint line




Endo-Model Modular Stems, uncemented

Modular Stems, uncemented, conical

MAT Tilastan – S

		Assembly length**							
		Tibia	Femur						
255	L	LT	LF1	LF2	LF3	LF4			
REF	mm	mm	mm	mm	mm	mm			
15-2952/01	50	72	89	92	96	99			
15-2952/02	80	102	119	122	126	129			
15-2952/03	95	117	134	137	141	144			
15-2952/04	120	142	159	162	166	169			
15-2952/05	135	157	174	177	181	184			
15-2952/06	160	182	199	202	206	209			
15-2952/07	200	222	239	242	246	249			
15-2952/08	240	262	279	282	286	289			
15-2952/09	280	302	319	322	326	329			






Modular Stems, uncemented, cylindrical

MAT Tilastan – S

				Assembly length**				
				Tibia				
	L	ØΑ	ØΚ	LT	LF1	LF2	LF3	LF4
REF	mm	mm	mm	mm	mm	mm	mm	mm
15-2951/01	60	10	16	82	99	102	102	109
15-2951/02	60	12	16	82	99	102	102	109
15-2951/03	60	14	16	82	99	102	102	109
15-2951/04	60	16	16	82	99	102	102	109
15-2951/05	60	18	18	82	99	102	102	109
15-2951/06	120	12	16	142	159	162	162	169
15-2951/07	120	14	16	142	159	162	162	169
15-2951/08	120	16	16	142	159	162	162	169
15-2951/09	120	18	18	142	159	162	162	169
15-2951/10	160	12	16	182	199	202	202	209
15-2951/11	160	14	16	182	199	202	202	209
15-2951/12	160	16	16	182	199	202	202	209
15-2951/13	160	18	18	182	199	202	202	209
15-2951/14	200	12	16	222	239	242	242	249
15-2951/15	200	14	16	222	239	242	242	249
15-2951/16	200	16	16	222	239	242	242	249
15-2951/17	200	18	18	222	239	242	242	249
15-2951/18	240	12	16	262	279	282	282	289
15-2951/19	240	14	16	262	279	282	282	289
15-2951/20	240	16	16	262	279	282	282	289
15-2951/21	240	18	18	262	279	282	282	289
15-2951/22	280	12	16	302	319	322	322	329
15-2951/23	280	14	16	302	319	322	322	329
15-2951/24	280	16	16	302	319	322	322	329
15-2951/25	280	18	16	302	319	322	322	329





Endo-Model Femoral Segments Tilastan,

for Endo-Model Modular Knee System – Rotating and Pure (Non-Rotating) Hinge Version



Femoral Segments, for Femoral Components: Size 1 and 2 MAT Tilastan – S

	For Femor	ral Components	:
REF	Size	Version	Width mm
Set: Size 1	(Height 20 mm)		\rightarrow
15-2971/00	x-small	right	55
15-2971/01	small	right	60
15-2971/02	medium	right	65
15-2971/03	large	right	75
15-2971/95	x-small	left	55
15-2971/10	small	left	60
15-2971/11	medium	left	65
15-2971/12	large	left	75

	For Femo			
REF	Size	Version	Width mm	REF
Set: Size 2	(Height 25 mm)			consisting of:
15-2972/00	x-small	right	55	15-2972/98 15-2972/99
15-2972/01	small	right	60	15-2972/04 15-2972/05
15-2972/02	medium	right	65	15-2972/06 15-2972/07
15-2972/03	large	right	75	15-2972/08 15-2972/09
15-2972/95	x-small	left	55	15-2972/96 15-2972/97
15-2972/10	small	left	60	15-2972/14 15-2972/15
15-2972/11	medium	left	65	15-2972/16 15-2972/17
15-2972/12	large	left	75	15-2972/18 15-2972/19

REF	Side
consisting of:	
15-2972/98	medial
15-2972/99	lateral
15-2972/04	medial
15-2972/05	lateral
15-2972/06	medial
15-2972/07	lateral
15-2972/08	medial
15-2972/09	lateral
15-2972/96	medial
15-2972/97	lateral
15-2972/14	medial
15-2972/15	lateral
15-2972/16	medial
15-2972/17	lateral
15-2972/18	medial
15-2972/19	lateral



Endo-Model Femoral Segments UHMWPE,

for Endo-Model Modular Knee System – Rotating and Pure (Non-Rotating) Hinge Version



Femoral Segments, for Femoral Components: Size 1 and 2

	For Femo	For Femoral Components:					
REF	Size	Version	Width mm	REF			
Set: Size 1	(Height 20 mm)		\rightarrow	consisting			
15-2965/01	x-small	right	55	15-2965/0 15-2965/0			
15-2961/01	small	right	60	15-2961/0 15-2961/0			
15-2961/02	medium	right	65	15-2961/0 15-2961/0			
15-2961/03	large	right	75	15-2961/0 15-2961/0			
15-2965/10	x-small	left	55	15-2965/1 15-2965/1			
15-2961/10	small	left	60	15-2961/1 15-2961/1			
15-2961/11	medium	left	65	15-2961/1 15-2961/1			
15-2961/12	large	left	75	15-2961/1 15-2961/1			


REF	Side
consisting of:	
15-2965/02	medial
15-2965/03	lateral
15-2961/04	medial
15-2961/05	lateral
15-2961/06	medial
15-2961/07	lateral
15-2961/08	medial
15-2961/09	lateral
15-2965/12	medial
15-2965/13	lateral
15-2961/14	medial
15-2961/15	lateral
15-2961/16	medial
15-2961/17	lateral
15-2961/18	medial
15-2961/19	lateral

	For Femo	ral Components			
REF	Size	Version	Width mm	REF	Side
Set: Size 2	2 (Height 25 mm)		\rightarrow	consisting of:	
15-2966/01	x-small	right	55	15-2966/02 15-2966/03	medial lateral
15-2962/01	small	right	60	15-2962/04 15-2962/05	medial lateral
15-2962/02	medium	right	65	15-2962/06 15-2962/07	medial lateral
15-2962/03	large	right	75	15-2962/08 15-2962/09	medial lateral
15-2966/10	x-small	left	55	15-2966/12 15-2966/13	medial lateral
15-2962/10	small	left	60	15-2962/14 15-2962/15	medial lateral
15-2962/11	medium	left	65	15-2962/16 15-2962/17	medial lateral
15-2962/12	large	left	75	15-2962/18 15-2962/19	medial lateral



Endo-Model Femoral Segments Tilastan,

for Endo-Model Modular Knee System - Rotating and Pure (Non-Rotating) Hinge Version



Femoral Segments, Size 3*

MAT Tilastan, UHMWPE

			For Femora	al Components:			
REF	Version	A mm	Size	Width mm	A mm	Version	REF
Size 3* (He	eight 50 mm)						
15-2973/01	right	60	small	60	60	left	15-2973/02
15-2973/03	right	65	medium	65	65	left	15-2973/04
15-2973/05	right	75	large	75	75	left	15-2973/06

* Only to be used in combination with longer stems (stem length above segments at least 180 mm). Not compatible with LINK MEGASYSTEM-C – Modular Joint Components Endo-Model with female taper.



Femoral Segments, Size 4*

MAT Tilastan

			For Femoral Components:				
REF	Version	A mm	Size	Width mm	A mm	Version	REF
Size 4* (He	ight 80 mm)						
15-2977/01	right	60	small	60	60	left	15-2977/02
15-2978/01	right	65	medium	65	65	left	15-2978/02
15-2979/01	right	75	large	75	75	left	15-2979/02

* Only to be used in combination with longer stems (stem length above segments at least 180 mm). Not compatible with LINK MEGASYSTEM-C – Modular Joint Components Endo-Model with female taper.



Endo-Model Femoral Segments Tilastan-S,

for Endo-Model Modular Knee System – Rotating and Pure (Non-Rotating) Hinge Version



Distal Femoral Segments (only to be used in combination with size 4) **MAT** Tilastan – S

REF	Height mm	for size
15-2977/10	10	small
15-2977/20	20	small
15-2977/40	40	small
15-2977/60	60	small
15-2977/80	80	small
15-2978/10	10	medium
15-2978/20	20	medium
15-2978/40	40	medium
15-2978/60	60	medium
15-2978/80	80	medium
15-2979/10	10	large
15-2979/20	20	large
15-2979/40	40	large
15-2979/60	60	large
15-2979/80	80	large

* Not compatible with LINK MEGASYSTEM-C – Modular Joint Components Endo-Model with female taper.





Endo-Model Femoral Segments UHMWPE,


for Endo-Model Modular Knee System - Rotating and Pure (Non-Rotating) Hinge Version



Femoral Segments, size 3*

			For Femoral Components:				
REF	Version	A mm	Size	Width mm	A mm	Version	REF
Size 3* (He	ight 50 mm)						
15-2967/01	right	55	x-small	55	55	left	15-2967/10
15-2963/01	right	60	small	60	60	left	15-2963/02
15-2963/03	right	65	medium	65	65	left	15-2963/04
15-2963/05	right	75	large	75	75	left	15-2963/06

* Only to be used in combination with longer stems (stem length above segments at least 180 mm). Not compatible with LINK MEGASYSTEM-C – Modular Joint Components Endo-Model with female taper.



Femoral Segments, size 4*

			For Femoral Components:				
REF	Version	A mm	Size	Width mm	A mm	Version	REF
Size 4* (He	ight 80 mm)						
15-2964/99	right	55	x-small	55	55	left	15-2964/00
15-2964/01	right	60	small	60	60	left	15-2964/02
15-2964/03	right	65	medium	65	65	left	15-2964/04
15-2964/05	right	75	large	75	75	left	15-2964/06

* Only to be used in combination with longer stems (stem length above segments at least 180 mm). Not compatible with LINK MEGASYSTEM-C – Modular Joint Components Endo-Model with female taper.



Endo-Model Femoral Segments UHMWPE,

for Endo-Model Modular Knee System – Rotating and Pure (Non-Rotating) Hinge Version



Femoral Segments

(only to be used in combination with size 4)

REF	Height mm	Size
15-2970/10	10	1
15-2970/20	20	2



Endo-Model Proximal Tibial Spacers UHMWPE,

for Endo-Model Modular Knee System – Rotating and Pure (Non-Rotating) Hinge Version



Proximal Tibial Spacers - straight -

MAT UHMWPE

REF	Size	Height REF Size mm
Set:		consisting of:
15-2516/70	x-small	15-2516/55x-small515-2516/60x-small1015-2516/65x-small15
15-2516/29	small	15-2516/05small515-2516/10small1015-2516/15small15
15-2517/29	medium	15-2517/05medium515-2517/10medium1015-2517/15medium15
15-2519/29	large	15-2519/05large515-2519/10large1015-2519/15large15

* Important Information: Proximal tibial spacers – straight – must not be combined with each other!

Proximal Tibial Spacers - anatomical -

MAT UHMWPE, EndoDur – S (CoCrM)

REF	Size	A Width mm	B Width mm
15-2516/24	x-small	55	40
15-2516/25	small	60	40
15-2517/26	medium	65	45
15-2519/27	large	75	55







Endo-Model Proximal Tibial Spacers Tilastan-S,

for Endo-Model Modular Knee System – Rotating and Pure (Non-Rotating) Hinge Version







Proximal Tibial Spacers, full right and left, incl. 2 hexagon socket countersunks with flat head screw 2.5 mm, for lateral and medial application, MAT Tilastan – S

REF	Size	H Height mm
15-2615/05	x-small	5
15-2615/10	x-small	10
15-2615/15	x-small	15
15-2616/05	small	5
15-2616/10	small	10
15-2616/15	small	15
15-2617/05	medium	5
15-2617/10	medium	10
15-2617/15	medium	15
15-2618/05	large	5
15-2618/10	large	10
15-2618/15	large	15



Proximal Tibial Spacers, half

incl. hexagon socket countersunks with flat head screw 2.5 mm, for lateral and medial application, MAT Tilastan – S

REF	Size	H Height mm
15-2990/11	x-small	5
15-2990/12	x-small	10
15-2990/13	x-small	15
15-2990/01	small	5
15-2990/04	small	10
15-2990/07	small	15
15-2990/02	medium	5
15-2990/05	medium	10
15-2990/08	medium	15
15-2990/03	large	5
15-2990/06	large	10
15-2990/09	large	15

Important Information:

Proximal tibial spacers of Tilastan must not be combined with each other!



Endo-Model Proximal Tibial Spacers Tilastan-S,

or Endo-Model Modular Knee System – Rotating and Pure (Non-Rotating) Hinge Version

Proximal Tibial Segments - anatomical -

MAT Tilastan – S

REF	Size	B Width mm	H Height mm
15-2981/01	x-small	55	50
15-2982/01	small	60	50
15-2983/01	medium	65	50
15-2984/01	large	75	50



Proximal Tibial Spacers

MAT Tilastan – S

REF	L Length mm	for Size
15-2981/10	10	x-small
15-2981/20	20	x-small
15-2981/40	40	x-small
15-2981/60	60	x-small
15-2982/10	10	small
15-2982/20	20	small
15-2982/40	40	small
15-2982/60	60	small
15-2983/10	10	medium
15-2983/20	20	medium
15-2983/40	40	medium
15-2983/60	60	medium
15-2984/10	10	large
15-2984/20	20	large
15-2984/40	40	large
15-2984/60	60	large







Centralizers, Patellar Components


Centralizers

MAT UHMWPE



Patellar Components centrical, circular MAT UHMWPE

REF	Size	Ø mm
15-2521/30	small	30
15-2521/35	medium	35
15-2521/40	large	40





Replacement Sets for Endo-Model Rotating Hinge and Pure (Non-Rotating) Hinge Version

Replacement Sets for Rotating Hinge Knee Prostheses, with anti-luxation device

MAT CoCrMo

REF	Side	Size
15-0027/10	right/left	x-small
15-0027/11	right/left	small
15-0027/12	right/left	medium
15-0027/13	right/left	large

Each package contains:

- complete coupling mechanism,
- bearing boxes,
- PE plateau and PE plateau anchoring screw.

Replacement Sets for LINK PorEx Rotating Hinge Knee Prostheses, with anti-luxation device MAT EndoDur – S (CoCrMo)/LINK PorEx*,UHMWPE

REF	Size
15-3027/10	x-small
15-3027/11	small
15-3027/12	medium
15-3027/13	large

Replacement Sets for LINK PorEx Rotating Hinge

Tibial Plateaus, with security screw MAT CoCrMo/LINK PorEx*, UHMWPE

UHMWPE/ CoCrMo REF	UHMWPE/ LINK PorEx* REF	Size
15-0027/17	15-0037/17	x-small
15-0027/14	15-0037/14	small
15-0027/15	15-0037/15	medium
15-0027/16	15-0037/14	large

Each package contains:

• PE plateau and PE plateau anchoring screw

Replacement Set for Endo-Model Rotational Bushing

MAT EndoDur – S (CoCrMo), UHMWPE

REF	Side	Size
15-1027/10	right/left	x-small
15-1027/11	right/left	small
15-1027/12	right/left	medium
15-1027/13	right/left	large

Each package contains: Rotational Bushing.

Replacement Sets for Pure (Non-Rotating) Hinge Knee Prostheses, with security screw

MAT EndoDur – S (CoCrMo), UHMWPE

Т

REF	Side	Size
15-0027/20	right	x-small
15-0027/21	right	small
15-0027/22	right	medium
15-0027/23	right	large
15-0027/30	left	x-small
15-0027/31	left	small
15-0027/32	left	medium
15-0027/33	left	large

Each package contains:

- Complete coupling mechanism,
- Bearing boxes,

• PE plateau and PE plateau anchoring screw.

* LINK PorEx: TiNbN = Titanium Niobium Nitride; surface modification (gold color). Available on request




MIRETO Instrument Set for Endo-Model Knee Prosthesis Systems

The **MIRETO Instrument Set** comprises seven instrument trays. It was developed with the focus on the following objectives:

- Instruments are securely held in their respective positions in the tray
- Ergonomic arrangement means the instruments are easy to remove, ensuring that surgery proceeds smoothly
- Optimally organized individual trays facilitate the job of the instrument nurse during the surgical procedure
- Trays are clearly marked with a picture of each instrument to ensure they can be quickly equipped with the correct instruments


On request, pictorial templates for equipping the instrument trays can be supplied. Cleaning instructions for all take-apart instruments are also available.



REF	MIRETO Instrument Set for Endo-Model Knee Joint Prostheses
15-6011/00	Case 1 – General Instruments
15-6012/00	Case 2 – Femoral Instruments (2 trays)
15-6013/00	Case 3 – Tibial Instruments
15-6014/00	Case 4 – Tapered Reamers conical & cylindrical
15-6015/00	Case 5 – Trial Prostheses (2 trays)
15-6016/00	Case 6 – Trial Prostheses cylindrical
15-6017/00	Case 7 – Tapered Reamers conical








15-6011/00 Case 1 – General Instruments

1	15-6001/00	Instrument Tray, Case 1, empty, 485 x 253 x 80 mm
2	64-1181/16	Hex Screwdriver, with metal handle, hex. 2 mm
3	10-5373/01	Hex Screwdriver, with metal handle, hex. 2.5 mm
4	15-8035/02	Inserter Positioner for PE plateaus
5	317-586	Inserter/Extraction Forceps for fixation pins Ø 3 mm, 210 mm
6	15-2537	Impactor Handle, small/medium, for femoral components
7	16-0115/01	Mallet, blow-back proof, incl. Spare Part Polyethylene
8	64-8008/02	Hex Screwdriver, hex 3.5 mm
9	317-658/01	Bone Awl with trocar point, 215 mm
10	15-2537/02	Impactor Handle, large, for femoral components
11	317-585/65 317-585/95 319-581/00	 Wire Pins, Ø 3 mm, 65 mm (4 pieces) Wire Pins, Ø 3 mm, 95 mm (4 pieces) optional Drill Pins, Ø 3 mm, 65/80 mm (4 pieces)
	319-582/00	Drill Pins, Ø 3 mm, 95/110 mm (4 pieces)
	15-6096/00	Alignment Device for modular stems with female taper (1 piece)
12	15-6053/00	T-Handle, with Hudson fitting
13	16-3287/00B	Adapter, for LINK power tool snap lock adapter, with Hudson fitting
14	16-3283/00 16-3284/00 16-3285/00	Adapter, optional with fitting: Hudson female/Jacobs male Hudson female/AO male Hudson female/Harris male
15	317-607/50	Cutting Template
16	15-6037/00	Drill, to open the femoral and tibial canal
17	16-3203/00	Impaction Plate for reamers
18	15-6098/00	Grooved Driver Tibia
19	15-6060/00	Stylus, for centralizer Ø 12 mm
20	15-6060/01	Stylus, for centralizer Ø 14 mm
21	15-6060/02	Stylus, for centralizer Ø 16 mm






1	15-6002/00	Instrument Tray, Case 2 (Tray 1), empty below, 485 x 253 x 80 mm			
		Basic Frame			
2	15-6030/01	Size S			
3	15-6030/02	Size M			
4	15-6030/03	Size L			
5	15-6111/00	Interlocking Pin for basic frame			
6	15-6110/00	Interlocking Spring for basic frame			
		Spacer, complete, 2 pieces			
7	15-6045/00	Height 2 mm			
8	15-6045/01	Height 10 mm			
9	15-6045/02	Height 15 mm			
10	15-6045/03	Height 20 mm			
11	15-6045/04	Height 25 mm			
12	15-6033/00	Alignment Rod, Ø 6 mm, 150 mm			
13	15-6040/01	Alignment Gauge, for sizes S/M/L			
14	15-6032/00	Drill for pin holes, Ø 6 mm			
15	15-6038/00	Depth Mill for taper coupling			
16	15-6049/00	Stylus, anterior			
17	15-6039/01	Saw Guide, for ventral rim, size M/L			
18	15-6039/00	Saw Guide, for ventral rim, size XS/S			
19	15-6042/00	Reamer, for box profile milling, Ø 12 mm, 74 mm			
20	15-6034/00	Slide-In Module, for ventral bone rim			
21	15-6046/00	Protective Cap, Ø 12 mm, 54 mm			
22	15-6046/01	Protective Cap, Ø 14 mm, 54 mm			
23	15-6031/00	Alignment Insert, right			
24	15-6031/01	Alignment Insert, left			

15-6012/00 Case 2 – Femoral Instruments (Tray 1)



15-6012/00 Case 2 – Femoral Instruments (Tray 2)



1	15-6002/00	Instrument Tray, Case 2 (Tray 2), empty above, 485 x 253 x 80 mm		
2	15-6036/00	Reamer for box pre-milling, Ø 24 mm		
3	15-6042/01	Reamer for box profile milling, Ø 12 mm, 100 mm		
		Basic Milling Guidances for condyle milling, complete		
4	15-6043/01	Size S		
5	15-6043/02	Size M		
6	15-6043/03	Size L		
7	15-6044/02	Reamer for condyle milling, Ø 26 mm		
8	15-6044/01	Inserting Forceps for condyle milling guidance		
9	15-6044/00	Condyle Milling Guidance for condyle milling		
10	15-6044/20	Screw for condyle milling guidance		
11	15-6044/30	Lever for condyle milling guidance		
12	15-6041/00	Slide-In Module for box profile milling (3 parts)		
13	15-6035/00	Slide-In Module for box pre-milling Endo-W		
14	15-6035/01	Slide-In Module for box pre-milling Endo-M Standard		



Г 17 16

1	15-6003/00	Instrument Tray, Case 3, empty, 485 x 253 x 80 mm			
2	16-3242/00	Alignment Rod Tibia			
3	15-6058/00	Tibial Saw Guide, 3 parts			
		Drill Templates Tibia			
4	15-6050/01	Size S			
5	15-6050/02	Size M			
6	15-6050/03	Size L			
7	16-3211/00	Guide Rod (2 pieces)			
8	15-6051/00	Alignment Gauge			
9	15-6054/01	Compressor, for Endo-Model-M, size XS/S			
10	15-6055/02	Compressor Extension, for Endo-Model, size M/L			
11	15-6054/02	Compressor, for Endo-Model-M, size M/L			
12	15-6056/02	Stem Compressor, size M/L			
13	15-6056/01	Stem Compressor, size S			
14	16-3197/00	T-Handle			
15	15-6052/00	Drill, Ø 20 mm			
16	16-3271/20	Drill Guide			
		Spacer, Tibial Alignment, for sizes S/M/L			
17	15-6059/00	Height 5 mm			
18	15-6059/01	Height 10 mm			
19	15-6059/02	Height 15 mm			
20	15-6059/03	Height 25 mm			
21	15-6062/00	Impaction Plate Tibia			
22	15-6057/00	Alignment Gauge, for tibia resection (3 parts)			

15-6013/00 Case 3 – Tibial Instruments

Instruments




1 2 3 4 5 1 1 2 3 4 5

1	15-6004/00	Instrument Tray, Case 4, empty, 485 x 253 x 80 mm			
		Tapered Rea	Tapered Reamers, cylindrical, with Hudson fitting B		
2	15-6048/00	Ø 12 mm			
3	15-6048/01	Ø 14 mm			
4	15-6048/02	Ø 16 mm			
5	15-6048/03	Ø 18 mm			
		Tapered Rea	mers, conical, with Hudson fitting B		
6	15-6047/01	Ø 16 mm	50 mm		
7	15-6047/02	Ø 16 mm	80 mm		
8	15-6047/03	Ø 16 mm	95 mm		
9	15-6047/04	Ø 16 mm	120 mm		
10	15-6047/05	Ø 16 mm	135 mm		
11	15-6047/06	Ø 16 mm	160 mm		
12	15-6047/07	Ø 16 mm	200 mm		
13	15-6047/08	Ø 16 mm	240 mm		
14	15-6047/09	Ø 16 mm	280 mm		

15-6014/00 Case 4 – Tapered Reamers conical & cylindrical






1	15-6005/00	Instrument Tray, Case 5 (Tray 1), empty below, 485 x 253 x 80 mm			
2 3 4	15-6065/01 15-6065/02 15-6065/03	Tibia Trial Prostheses, intracondylar Size S Size M Size L			
5 6 7 8 9 10	15-6067/01 15-6067/02 15-6067/03 15-6068/01 15-6068/02 15-6068/03	Femur Trial Prostheses, intracondylarSize SleftSize MleftSize LleftSize SrightSize LrightSize Lright			
11 12 13 14	15-6088/01 15-6088/02 15-6093/01 15-6093/02	Femoral Trial Segments, universal Height 20 mm right Height 20 mm left Height 25 mm right Height 25 mm left			
15 16 17	15-6061/00 15-6094/00 15-6070/00	Extraction Instrument for trial prostheses (2 parts) Trial Support Ring, height 10 mm, Ø 28 mm Femoral Trial Adapter, for modular stems M10/M12			
18 19 20	15-6066/01 15-6066/00 15-6066/20	Trial Connection Component for hinged knee versions sizes S/M/L Trial Connection Component for rotational knee versions sizes S/M/L Trial Axis, for sizes S/M/L			
21 22 23 24 25	15-6078/03 15-6079/03 15-6080/03 15-6081/03 15-6082/03	Trial Segments Tibia Height 5 mm Sizes L right Height 5 mm Sizes L left Height 10 mm Sizes L right Height 10 mm Sizes L left Height 15 mm Sizes L left			
26 27 28 29 30 31 32	15-6083/03 15-6078/02 15-6079/02 15-6080/02 15-6081/02 15-6082/02 15-6083/02	Height 15 mmSizes LleftHeight 5 mmSize MrightHeight 5 mmSize MleftHeight 10 mmSize MrightHeight 10 mmSize MleftHeight 15 mmSize MrightHeight 15 mmSize MleftHeight 15 mmSize Mleft			
33 34 35 36 37 38	15-6078/01 15-6079/01 15-6080/01 15-6081/01 15-6082/01 15-6083/01	Height 15 mmSize MleftHeight 5 mmSize SrightHeight 5 mmSize SleftHeight 10 mmSize SrightHeight 10 mmSize SleftHeight 15 mmSize SleftHeight 15 mmSize SleftHeight 15 mmSize Sleft			

15-6015/00 Case 5 – Trial Prostheses (Tray 1)




15-6015/00 Case 5 – Trial Prostheses (Tray 2)



1	15-6005/00	Instrument Tray, Case 5 (Tray 2), empty, 485 x 253 x 80 mm			
		Tibial Trial Stems, for Endo-Model prosthesis stems			
2	15-6063/01	Size S			
3	15-6063/02	Size M			
4	15-6063/03	Size L			
		Femoral Trial Stems, for Endo-Model prosthesis stems			
5	15-6064/01	Size S			
6	15-6064/02	Size M			
7	15-6064/03	Size L			
		Trial Stems, conical, for femoral and tibial components, cemented/uncemented			
8	15-6071/01	Length 50 mm			
9	15-6071/02	Length 80 mm			
10	15-6071/03	Length 95 mm			
11	15-6071/04	Length 120 mm			
12	15-6071/05	Length 135 mm			
13	15-6071/06	Length 160 mm			
14	15-6071/07	Length 200 mm			
15	15-6071/08	Length 240 mm			
16	15-6071/09	Length 280 mm			





15-6016/00 Case 6 – Cylindrical Trial Stems

1	15-6006/00	Instrument Tra	y, Case 6, empty, 485 x 253 x 80 mm	
		Trial Stems, cylindrical, uncemented		
2	15-6073/01	Ø 12 mm	Length 120 mm	
3	15-6073/02	Ø 14 mm	Length 120 mm	
4	15-6073/03	Ø 16 mm	Length 120 mm	
5	15-6073/04	Ø 18 mm	Length 120 mm	
6	15-6074/01	Ø 12 mm	Length 160 mm	
7	15-6074/02	Ø 14 mm	Length 160 mm	
8	15-6074/03	Ø 16 mm	Length 160 mm	
9	15-6074/04	Ø 18 mm	Length 160 mm	
10	15-6075/01	Ø 12 mm	Length 200 mm	
11	15-6075/02	Ø 14 mm	Length 200 mm	
12	15-6075/03	Ø 16 mm	Length 200 mm	
13	15-6075/04	Ø 18 mm	Length 200 mm	
14	15-6076/01	Ø 12 mm	Length 240 mm	
15	15-6076/02	Ø 14 mm	Length 240 mm	
16	15-6076/03	Ø 16 mm	Length 240 mm	
17	15-6076/04	Ø 18 mm	Length 240 mm	
18	15-6077/01	Ø 12 mm	Length 280 mm	
19	15-6077/02	Ø 14 mm	Length 280 mm	
20	15-6077/03	Ø 16 mm	Length 280 mm	
21	15-6077/04	Ø 18 mm	Length 280 mm	
22	15-6072/00	Ø 10 mm	Length 60 mm	
23	15-6072/01	Ø 12 mm	Length 60 mm	
24	15-6072/02	Ø 14 mm	Length 60 mm	
25	15-6072/03	Ø 16 mm	Length 60 mm	
26	15-6072/04	Ø 18 mm	Length 60 mm	



1 8 10 12 14	13 11 9
6 4 2	3 5 7

15-6017/00 Case 7 – Conical Tapered Reamers

1	15-6007/00	Instrument Tray, Case 7, empty, 485 x 253 x 80 mm			
		Tapered Ream	Tapered Reamers: conical, with fitting B: Hudson		
2	16-5130/12	Ø 12 mm	Stem length 130 mm		
3	16-5130/13	Ø 13 mm	Stem length 130 mm		
4	16-5130/14	Ø 14 mm	Stem length 130 mm		
5	16-5130/15	Ø 15 mm	Stem length 130 mm		
6	16-5130/16	Ø 16 mm	Stem length 130 mm		
7	16-5130/17	Ø 17 mm	Stem length 130 mm		
8	16-5130/18	Ø 18 mm	Stem length 130 mm		
9	16-5130/19	Ø 19 mm	Stem length 130 mm		
10	16-5130/20	Ø 20 mm	Stem length 130 mm		
11	16-5130/21	Ø 21 mm	Stem length 130 mm		
12	16-5130/22	Ø 22 mm	Stem length 130 mm		
13	16-5130/23	Ø 23 mm	Stem length 130 mm		
14	16-5130/24	Ø 24 mm	Stem length 130 mm		



Additional Instrument Set

for Endo-Model Rotating Hinge Knee and Endo-Model Modular Knee Prostheses



REF	for Endo-Model Rotating Hinge Knee and Modular Knee Prostheses	
15-2529/90	Set, complete, in 1 small container K1, on 1 tray with storage racks	
05-1000/01	Small Container K1, only, 460 x 190 x 92 mm	
15-2529/91	Tray, empty, 405 x 165 x 50 mm	
64-8008/02	Hex Screwdriver with metal handle, hex 3.5 mm, 250 mm	
15-2544	 Separate Rod for removal of the rotating bushing version V02, Ø M5, 210 mm Hex Screwdriver with metal handle, hex 2.5 mm, 180 mm Torque Wrench, hex 2.5 mm, 205 mm 	
10-5373/01		
15-2545		



Tibial Preparation



01

Mark the entry point and open the tibial canal with the drill (15-6037/00) at the point of attachment of the anterior cruciate ligament. The drill can either be used manually with the T-handle (15-6053/00) or driven mechanically.

02

Assemble the tapered reamer with the T-handle (15-6053/00). The tibial stop plate (15-6062/00) is clicked into the designated groove on the tapered reamer stem. When using uncemented modular stems: Ream with increasing diameters until the awl has achieved contact with the cortical bone over a continuous distance of approx. 50 mm. The implant must correspond with the tapered reamer last used in terms of diameter and length.

In the case of cemented modular stems, the tapered reamer should be at least 2 mm larger than the planned stem diameter.

Important information:

The position of the stop plate represents the lower edge of the tibial metal tray. The tapered reamers must not be used with a drive.







Once the tapered reamer is securely positioned by means of contact with the cortical bone, the T-handle and the stop plate are both removed. The tapered reamer remains in the medullary space. The tibial resection alignment instrument (15-6057/00) is attached, in a pre-assembled state (scale setting "0"), onto the tapered reamer. It must be ensured here that the stylus pin, which can be positioned medially or laterally, is guided in the designated groove.

04

The saw guide (15-6058/00) is attached to the alignment guide and placed on the ventral tibia. Ideally, the saw guide is pre-assembled in such a way that tibial resection can be performed from medial.







By turning the setting screw on the alignment guide, the desired resection height can be set and read off the scale. The cutting template (317-607/50) can be used to check the resection height. For first-time surgery, the resection height should not be greater than 10 mm. For revision surgery, as far as possible very little or no bone should be additionally resected. Resection is performed without a dorsal slope at a 90° angle to the tibial axis (default setting due to instrument design).

06

After setting the desired resection height, the saw guide is fixed in place using two fixation pins (317-585/65 or /95) or drill pins (319-581/00 or 319-582/00). For this, the pins are initially inserted into the holes on the line marked "0" in the medial holes groups (results in parallel alignment of the pins). After removing the tapered reamer, a third fixation pin or drill pin is inserted at an angle in one of the neighboring hole groups to secure the saw guide in place.





After removing the tapered reamer, the bone can be resected. To achieve an optimal cutting result, saw blades with a width of 1.24 - 1.27 mm are recommended. The fixation pin inserted at an angle is then removed and the saw guide is pulled off toward the front. The two parallel fixation or drill pins can remain in place. These allow the saw guide to be reattached at the +2 or +4 mm level in order to repeat resection later on.

80

The tapered reamer used last is again inserted in the medullary canal with the stop plate attached. The final implant size is determined by applying a drill template (15-6050/01, /02, /03) which corresponds exactly in size with the respective implant. Since the Endo-Model implants are also anchored by the diaphysis, full cortical bone support of the prosthesis is not necessary. The implant must not project over the bone margin.

The alignment guide (15-6051/00) is guided over the stem of the tapered reamer and connected to the cylindrical elevations of the drill template.

For revision surgery, the level of the tibial preparation can be raised to the level of the planned and to be reconstructed joint line. For this, spacers (15-6059/00, /01, /02, /03) of the relevant implant heights can be pushed on the drill template.







For rotational alignment, the alignment rod (16-3242/00) can be pushed from ventral onto the drill template. The ventral alignment rod must be aligned in the area between the middle of the tibial tuberosity and its medial margin.

10

After rotational alignment, the drill template is fixed with four fixation pins or drill pins. Mark the rotary alignment on the ventral cortical bone of the tibial head.







Remove the alignment guide, alignment rod and tapered reamer.

12

Attach the drill guide (16-3271/20) to the drill template. The lower collar must sit flush. The collarless area points toward dorsal.







Insert the drill (16-6052/00) into the proximal tibia as far as it will go.

14

Screw the guide rods (16-3211/00) into the dorsal thread holes of the drill template.



MIRETO Surgical Technique





- Attach the compressor addition (15-6055/02) onto the compressor (15-6054/02).
- Screw the stem compressor (15-6056/02) onto the compressor and attach the handle.

16

Insert the compressor via the guide rods as far as possible on the drill template. Finally remove all instruments.







Assemble the tibial trial. If necessary, push in trial washers from medial and/or lateral.

18 Insert the trial prosthesis.





Femoral Preparation



Mark the entry point and open the femoral canal with the drill (15-6037/00). The drill can either be used manually with the T-handle (15-6053/00) or driven mechanically.

02

Assemble the tapered reamer with the T-handle (15-6053/00). The femoral stop plate (16-3203/00) is clicked into the designated groove on the tapered reamer stem. When using uncemented modular stems: Ream with increasing diameters until the awl has achieved contact with the cortical bone over a continuous distance of approx. 50 mm. The implant must correspond with the tapered reamer last used in terms of diameter and length.

In the case of cemented modular stems, the tapered reamer should be at least 2 mm larger than the planned stem diameter.

Once the tapered reamer is securely positioned by means of contact with the cortical bone, the T-handle and the stop plate are both removed. The tapered reamer remains in the medullary canal.

Important information:

The position of the stop plate represents the level of the joint line. The tapered reamers must not be used with a drive.








Assemble the base frame (15-6030/01, /02, /03) in line with the size of the tibial preparation:

- Insert the locking pin (15-6111/00) into the groove on the relevant side L (left) or R (right).
- Insert the locking spring (15-6110/00) into the outer drill hole opposite the locking pin.
- The alignment insert (15-6031/00, /01) is then attached to the front of the base frame and secured in place by closing the locking pin.



04

Attach the base frame onto the stem of the tapered reamer and guide the instrument on the joint surface.






With revision surgery or extreme varus or valgus deformities, different or same-sized spacers (15-6045/00 to/04) can be inserted laterally into the base frame medially, laterally or on both sides for a more stable fixation and distalization of the joint line.

06

Check the emergence point of the angled anterior saw cut. Insert the alignment rod (15-6033/00) as far as possible into the alignment insert. Position the anterior stylus (16-6049/00) in the designated groove in the alignment rod and push up to the ventral bone. If the marking line of the selected prosthesis size is around the middle of the alignment rod, the saw blade will emerge precisely at the contact point of the ventral cortical bone. If the marking is significantly below the lower third or below the alignment rod, the position of the tapered reamer must be corrected to ventral and the emergence point checked again.








Check the incision for femoral segments (optional) and the anterior and posterior saw cuts using the cutting template (317-607/50).

80

Insert the alignment guide (15-6040/01) into the slot of both cylindrical drill holes in the base frame as far as possible. Check and determine the rotation/ position of the prosthesis axis with respect to the epicondyles. The position of the tips of the alignment guide represents the distance of the prosthesis axis in relationship to the position of the joint line.









Fix the base frame by inserting four fixation pins (317-585/65 or /95) or drill pins (319-581/00 or 319-582/00) in the medial and lateral drill holes.

10

Remove the alignment insert and the tapered reamer. Attach and secure the slide-in unit for the ventral bone margin (15-6034/00) in the T-groove of the base frame from above. Insert the short reamer (15-6042/00) three times as far as it will go.






Pre-drill two recesses for the medial and lateral fixation pegs for the Endo-Model–M or –W knee prostheses using the drill for making drill holes for pins (15-6032/00). This step is not necessary for the Endo-Model Standard prosthesis.




12

Remove the previous insert. Attach and secure the saw guide for the ventral margin (15-6039/00, /01) of the same size in relation to the base frame. Make the anterior saw cut. To achieve an optimal cutting result, saw blades with a width of 1.24 - 1.27 mm are recommended.








Make the posterior saw cuts through the two medial and lateral saw slots in the base frame. The saw guide is to remain in the base frame so as to stop the saw blade from slipping out of the guide and inwards.

Endo Std/M

14

Remove the saw guide. Insert and secure the slide-in unit for box preliminary milling (15-6035/00, /01). The choice of slide-in unit depends on the Endo-Model prosthesis used (Standard-M or Standard-W). Insert the reamer twice as far as it will go for preliminary milling (15-6036/00).





Pre-drill the cylindrical part of the tapered coupling when using the Endo-Model–W knee prosthesis. Attach and secure the alignment insert (15-6031/00, /01). Insert the drill (15-6038/00) in the alignment insert with the narrowest part of the drill shank. With the Endo-Model Standard and Endo-Model Standard–M prosthesis this step is not necessary.



16

Assemble the slide-in unit for box preliminary milling (15-6041/00) as described in fig. **A** and **B**. The black marking lines must be parallel when assembled.







For preliminary milling of the box: Attach and secure the slide-in unit from above in the T-groove of the base frame. Starting with the short reamer (15-6042/00) – inserted as far as it will go – the remaining intracondylar bone is milled in a clockwise direction. The individual steps are then repeated with the long reamer (15-6042/01). The bone must be prepared carefully to ensure that the next instrument fits perfectly when inserted.



18

Pre-assemble the base reamer guide to mill the condylar (15-6043/01, /02, /03). Insert the U-shaped lever and tighten the screw using the 3.5 mm hex socket screwdriver (64-8008/02).








Connect the base reamer guide with the inserting forceps (15-6044/01) to mill the condylar.

20


Intracondylar insertion of the base reamer guide to mill the condylar. The inserting forceps are guided through the recesses in the base frame and, when finally positioned, their front ends must rest against the frame (see arrow).











Insert the reamer for condylar milling (15-6044/02) into the left-hand side of the reamer guide. Once started, the reamer is guided as far as it will go and moved from dorsal to ventral. It may be necessary to repeat this step several times until no more cartilage is removed. The reamer is then reset and the individual steps are repeated on the right-hand side. Here work is from ventral to dorsal.

25

Remove the instruments for condylar milling.





Trial Reduction



26

Assemble and insert the trial femoral component.

27

Insert the trial axis (15-6066/20) in the trial connecting component for the rotational and hinged versions (15-6066/00, /01) and push onto the trial tibial component as far as possible.









Connect the two joint pieces by inserting the dorsal recess of the trial axis into the axis of the femoral component and then pushing the tibia up. Test the prosthesis. Deficits in the flexion and extension gap are compensated for with trial femoral segments and/or trial tibial washers. The components are separated by proceeding in reverse order.

29

Remove the trial femoral and tibial components with the extraction instrument for trial prostheses (15-6061/00).





Implantation



30

Before inserting the cemented stems, check the depth and suitable diameter of the centering stars (12, 14 or 16 mm) with the stylus for centering stars (15-6060/00, /01, /02).

31

After assembling the suitable centering star (with cemented stems), the femoral component is inserted using the impactor for femoral components (15-2537, 15-2537/02). Then any excess cement and the plastic connector (prevents cement entering the joint mechanism) are removed. The connector is to be pulled toward dorsal. With full cementation, as a rule we recommend using at least 80 g of bone cement.






After assembling the appropriate centering star (for cemented stems), the tibial component is positioned using the tibial driver (15-6098/00). When aligning the prosthesis, the previously marked rotation setting must be taken into account.

Important:

The tibial component may only be implanted without a polyethylene plateau with a trial screw screwed in as far as possible. To remove the polyethylene plateau, the trial screw is undone with the 3.5 mm hex socket screwdriver (64-8008/02) and the polyethylene plateau is removed using the inserter (15-8035/02). After that, the trial screw is screwed in as far as it will go. This is the only way to prevent cement entering the drill hole in the metal plateau. Any excess cement is removed. For complete cementation, as a rule we recommend using at least 40 g of bone cement.

33

The polyethylene plateau also has to be removed from the tibial metal tray when the components are fitted together. The lower part pin is inserted into the femoral component with the knee flexed.





To insert the polyethylene tibial plateau, the femoral component is lifted slightly. The tibial plateau is then inserted from ventral between the upper and lower part of the prosthesis. It must be ensured that the chamber of the plastic plateau grips over the flange of the femoral component and that the dovetail-shaped incision on the bottom of the UHMWPE plateau fits into the peripheral groove on the metal tibial support.



35

The PE plateau is secured into place on the tibial metal tray using the 3.5 mm hex socket screwdriver (64-8008/02) and the self-locking fixation screw.

Important:

A self-locking fixation screw may only be used for the final assembly of the plateau. If the fixation screw is undone, the screw lock in the PE plateau is destroyed and a new plateau must be inserted.

The implanted rotating hinge knee prosthesis should enable flexion of up to 90° providing that the soft tissue allows this. When extended, slight elastic extension restriction of around 5° is ideal. This serves to ensure the reliable closure of both prosthesis components.



Literature



E. Engelbrecht, A. Siegel, J. Röttger, and Prof. H. W. Buchholz* Statistics of Total Knee Replacement: Partial and Total Knee Replacement, Design St. Georg Journal of Clinical Orthopaedics, 1976, No. 120, pp 54-64 (K3)

E. Engelbrecht, E. Nieder, E. Strickle, A. Keller Intrakondyläre Kniegelenkendoprothese mit Rotationsmöglichkeit – ENDO-MODELL® *CHIRURG 52: 368-375 (1981)* (K1)

R. Dederich und L. Wolf Kniegelenkprothesen-Nachuntersuchungsergebnisse Unfallheilkunde (1982) 85:359-368 (K2)

J. Röttger, K. Heinert Die Knieendoprothesensysteme (Schlitten- und Scharnierprinzip). Beobachtungen und Ergebnisse nach 10 Jahren Erfahrung mit über 3700 Operationen. Z. Orthop. 122(1984) 818-826 (K17)

E. Nieder, E. Engelbrecht, A. Keller **Totale intrakondyläre Scharniergelenkendoprothese mit Rotationsmöglichkeit – Endo-Modell**[®] *Sonderdruck aus Heft 5: Orthopädische Praxis*, 1987, 23. Jahresgang, Seite 402-412 **(K34)**

K. Heinert, E. Engelbrecht Total Knee Replacement - Experience with a Surface and Total Knee Replacement: Further Development of the Model St. Georg[®]. 2400 Sledges and Hinges

Proceedings of the International Symposium on Total Knee Replacement, May 19-20, 1987, Nagoya, Japan Springer Verlag:, Berlin Heidelberg, New York Tokyo (1987), pp 257-273 **(K53)**

E. Engelbrecht, M.D. **The Tibial Rotating Knee Prosthesis "Endo" Model: Surg. Technique** *The Journal of Orthopaedic Surgical Techniques, Volume 3, Number 2, 1987* **(K36)**

K. Heinert, E. Engelbrecht Langzeitvergleich der Knie-Endoprothesensysteme St. Georg® 10-Jahres-Überlebensraten von 2236 Schlittenund Scharnier-Endoprothesen Der Chirurg (1988) 59:755-762 (K38)

F. Madsen, P. Kjarsgaard-Andersen, M. Juhl, O. Sneppen A Custom-Made Prosthesis for the Treatment of Supracondylar Femoral Fractures after Total Knee Arthroplasty: Report of Four Cases Journal of Orthopaedic Trauma, Vol. 3, No. 4, pp. 333-337,1989 (K42)

E. Nieder

Schlittenprothese, Rotationsknie und Scharnierprothese Modell St. Georg[®] und Endo-Modell[®]. Differentialtherapie in der primären Kniegelenkalloarthroplastik *Orthopäde (1991) 20:170-180* (K45)

G. von Förster, D. Klüber und U. Käbler Mittel- bis langfristige Ergebnisse nach Behandlung von 118 peripro-thetischen Infektionen nach Kniegelenkersatz durch einzeitige Austauschoperationen Orthopäde(1991) 20: 244-252 (K46)

Adolph V. Lombardi, Jr, Thomas H. Mallory, Robert W. Eberle, and Joanne B. Adams Results of Revision Total Knee Arthroplasty Using Constrained Prostheses Seminars in Arthroplasty, Vol 7, No. 4 (October), 1996: pp 349-355

E. Engelbrecht, E. Nieder, D. Klüber

Reconstruction of the Knee - Ten to Twenty Years of Knee Arthroplasty at the Endo-Klinik: A Report on the Long-term Follow-up of the St. Georg[®] Hinge and the Medium-term Follow-up of the Rotating Knee Endo-Model[®] Springer Verlag: Tokyo, Berlin, Heidelberg, New York (1997) (K57)

E. Nieder **Revisionsalloarthroplastik des Kniegelenks** Sonderausgabe aus: Orthopädische Operationslehre, Band II/1: Becken und untere Extremität Herausgegeben von R. Bauer, F. Kerschbaumer und S. Poisel F. Alt, U. Sonnekalb, N. Walker Unikondyläre Schlittenprothese versus scharniergeführte Totalendoprothesen des Kniegelenkes Orthopädische Praxis 1/98, 34. Jahresgang, Seite 20-24, 1998 (K61)

A. V. Lombardi, T. H. Mallory, R. E. Eberle, J. B. Adams **Rotating Hinge Prosthesis in Revision Total Knee Arthroplasty: Indications and Results** *A Reprint from Surgical Technology International VI, 1998* (K55)

E. Nieder, G.W. Baars, A. Keller **Totaler Tibia-Ersatz Endo-Modell**[®] *Orthopädie Aktuell: Nr. 5/1998, LINK News* (K60)

S. Schill, H. Thabe Die periprothetische Knieinfektion – Therapiekonzept, Wertigkeit und mittelfristige Ergebnisse Aktuelle Rheumatologie, Heft 5, 24. Jahrgang, 1999, pp 153-160 (K70)

G.W. Baars Knieendoprothetik: Das optimale Implantat für jeweilige Indikation finden Orthopäde 2000 (Suppl1) 29: S1-2

M. Zinck, R, Sellkau Rotationsknieprothese Endo-Modell®- Geführter Oberflächenersatz mit Sti(e)I Orthopäde 2000 (Suppl1) 29: S 38-42

M. Crowa, E. Cenna, C. Olivero Rotating knee prosthesis – Surface or hinge replacement? *Orthopäde 2000 (Suppl1) 29: S 43-44*

J-N. Argenson. J M. Aubaniac **Total Knee arthroplasty in femorotibial instability** *Orthopäde 2000.29:S 45-47, Springer Verlag 2000* **(K72)**

M. von Knoch, R. Brocks, C. Siegmüller, G. Ribaric, L. Leupolt, G. von Förster Knieflexion nach Rotationsknieendoprothese Z. Orthop 2000; 138: 66-68 (K71)

R.E. Windsor, K. Steinbrink Controversies in Total Knee Replacement Two-stage exchange is the optimal treatment for an infected total knee replacement Oxford University Press 2001 (K78)

Thomas Nau, MD, E. Pflegerl, MD, J. Erhart, MD, and V. Vecsei, MD Primary Total Knee Arthroplasty for Periarticular Fractures The Journal of Arthroplasty, Vol 18, No 8, 2003 (K82)

G. Petrou, H. Petrou, C. Tilkeridis, T. Stavrakis,
T. Kapetsis, N. Kremmidas, M. Gavras
Medium-term results with a primary cemented rotating-hinge total knee replacement
A 7-TO 15-YEAR FOLLOW-UP
J Bone Joint Surg (Br), 2004; 86-B :813-17 (K84)

M.R. Utting, J.H. Newman Customised hinged knee replacement as a salvage procedure for failed total knee arthroplasty *The Knee 11 (2004) 475-479* (K86)

Nayana Joshi, Antonio Navarro-Quilis Is There a Place for Rotating-Hinge Arthroplasty in Knee Revision Surgery for Aseptic Loosening? The Journal of Arthroplasty 2008; 23(8):1204-1210 (K94)

M. Napp, M. Frank, M. Witt Pathologische Fraktur des distalen Femurs bei Knie-TEP Der Orthopäde, Band 38, Heft 10, Oktober 2009 (K96)

Dae Kyung Bae, Sang Jun Song, Kyoung Ho Yoon, Jung Ho Noh Long-Term Outcome of Total Knee Arthroplasty in Charocot Joint: A 10- to 22- Year Follow-Up The Journal of Arthroplasty 2009; 24(8):1152-1156 (K98)







Please note the following regarding the use of our implants:

1. Choosing the right implant is very important.

The size and shape of the human bone determines the size and shape of the implant and also limits the load capacity. Implants are not designed to withstand unlimited physical stress. Demands should not exceed normal functional loads.

2. Correct handling of the implant is very important.

Under no circumstances should the shape of a finished implant be altered, as this shortens its life span. Our implants must not be combined with implants from other manufacturers. The instruments indicated in the Surgical Technique must be used to ensure safe implantation of the components.

3. Implants must not be reused.

Implants are supplied sterile and are intended for single use only. Used implants must not be used again.

4. After-treatment is also very important.

The patient must be informed of the limitations of the implant. The load capacity of an implant cannot compare with that of healthy bone!

- 5. Unless otherwise indicated, implants are supplied in sterile packaging.
- Note the following conditions for storage of packaged implants:
- Avoid extreme or sudden changes in temperature.
- Sterile implants in their original, intact protective packaging may be stored in permanent buildings up until the "Use by" date indicated on the packaging.
- They must not be exposed to frost, dampness or direct sunlight, or mechanical damage.
- Implants may be stored in their original packaging for up to 5 years after the date of manufacture. The "Use by" date is indicated on the product label.
- Do not use an implant if the packaging is damaged.

6. Traceability is important.

Please use the documentation stickers provided to ensure traceability.

7. Further information on the material composition is available on request from the manufacturer.

Follow the instructions for use!

Waldemar Link GmbH & Co. KG, Hamburg

All content in this catalog, including text, pictures and data, is protected by law. Every instance of use, whether in part or in whole and which is not permitted by law, is subject to our prior consent. In particular, this applies to the reproduction, editing, translation, publishing, saving, processing, or passing on of content stored in databases or other electronic media and systems, in any manner or form. The information in the catalogs is solely intended to describe the products and does not constitute a guarantee.

The Surgical Technique described has been written to the best of our knowledge and belief, but it does not relieve the surgeon of his/her responsibility to duly consider the particularities of each individual case.

Products shown in this document may not be available in your country. The product availability is subject to the approval and/or registration regulations of the respective country. Please contact Waldemar Link GmbH & Co. KG if you have questions about the availability of LINK products in your country.

Waldemar Link GmbH & Co. KG and/or other corporate affiliated entities own, use or have applied for the following trademarks in many jurisdictions: LINK, BiMobile, SP II, Modell Lubinus, E-Dur, EndoDur, T.O.P. II, BetaCup, CombiCup PF, CombiCup SC, CombiCup R, MobileLink, C.F.P., LCU, SP-CL, LCP, MIT-H, Endo-Model, Endo-Model SL, MP, MEGASYSTEM-C, GEMINI SL, SPAR-K, LCK, Link OptiStem, HX, TiCaP, X-LINKed, PorAg, LINK PorEx, BiPorEx, PorEx-Z, TrabecuLink, Tilastan, customLINK, RescueSleeve, Stactip, VACUCAST.

Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners.

This publication is intended for distribution in the U.S. only.

Waldemar Link GmbH & Co. KG

Barkhausenweg 10 · 22339 Hamburg · Germany Phone +49 40 53995-0 · info@link-ortho.com www.link-ortho.com

