




LinkSymphoKnee) Built on TRUST



"For my patients, the LinkSymphoKnee has solved many issues. The sizing fits all morphologic types without compromise. The system accommodates to the surgeon's preference rather than the surgeon having to accommodate to the system. The instrumentation is efficient and streamlined, ideal for both the hospital and the ASC."

- Dr. Matt Austin MD Rothman Orthopaedic Institute Philadelphia, PA





PS/PS+ Jump Height = 17mm

PS Constraint Free to Rotate

PS+ Constraint Constrained 3° IE and VV

Anterior Recess For Patella Ligament Wrapping

E-Dur Vitamin E Infused Polyethylene¹

5° Sloped Post Increases Stability

Posterior Chamfer Designed to Reduce Impingement

0° Posterior Slope

Posterior and Central
Dovetail Locking
Feature

